



AIRSSENTIAL HOME
Health Care Solutions

LifeTemp® Touch Forehead Thermometer



Model AI-TOU

Instruction Manual

Read the instruction manual carefully before using this device. Keep this manual in a safe place for future reference.

WARNING: Keep this device away from children and the infirm, to avoid accidents. Failure to follow the instructions may cause inaccurate readings or damage the thermometer.

Introduction

Congratulations on the purchase of your new LifeTemp® Touch Infrared Forehead Thermometer. The thermometer is simple to use, utilises German technology and sensors to ensure high accuracy, is rechargeable and has blue backlighting for use in dark environments.

What Forehead Temperature is a Fever?

An exact relationship between oral, ear, armpit, rectal and forehead temperature measurements has not yet been determined by medical research. By consensus a forehead temperature of 38.0° C (100.4 °F) or higher constitutes a fever.

Package Contents

This carton contains:
1 x Touch Forehead Thermometer
1 x Micro USB/USB charging cable
1 x User Instruction Manual

The LifeTemp Touch Forehead Thermometer
Thermometer Parts Description: The parts of this thermometer are identified in Figure 1.




FIGURE 1.

PARTS OF THE LIFETEMP TOUCH FOREHEAD THERMOMETER

LCD Display

The LifeTemp Touch LCD display is identified in Figure 2.




FIGURE 2. THE LIFETEMP TOUCH LCD DISPLAY.

Instructions for Use

- Turn the thermometer on by pressing On/Off switch (see Fig 3.1). The LCD screen will illuminate with a blue backlight.
- A beep sound will be heard and the backlight will turn off. The screen will then flash ---°C. The unit is now ready for use. Place the thermometer's probe on the forehead (see Fig. 3.2). After approximately 1 second the thermometer will beep and display the temperature with a blue backlight.
If a consecutive reading is required then please wait for the backlight to turn off and the °C to flash. Repeat the above procedure as often as required.
If a reading above 38°C is recorded then 6 consecutive beeps will be heard.
- Press the On/Off switch to turn the thermometer off (see Fig. 3.3). Or the unit will switch off automatically after approximately 60 seconds.



FIGURE 3.

Recalling Readings from Memory

- To enter memory mode, ensure the thermometer is Off.
- Touch and release thermometer's probe with your finger.
- The last reading will now be displayed as 01 and followed with an illuminated display. Touch and release the probe again to display the next reading in a sequence from last reading to first reading (see Fig.4).
- The thermometer will turn-off after 12 seconds of idling.
- The thermometer can recall 12 temperature readings.



FIGURE 4. RECALLING READINGS FROM MEMORY

Deleting Temperature Readings Stored in Memory

- Ensure the thermometer is Off.
- Cover the thermometer's probe for approximately 5 seconds. The LCD will display "dEL".
- Uncover and press the thermometer's probe again to delete all readings in memory (see Fig. 5).
- The thermometer will switch-off after 5 seconds of idling.



FIGURE 5. DELETING READINGS STORED IN MEMORY



Changing Temperature Scales

1. Ensure the thermometer is Off.
2. Cover the thermometer's probe for approximately 10 seconds to enter the temperature mode. The current temperature unit will flash.
3. Uncover and press the thermometer's probe again to toggle and select between °C / °F (see Fig. 6).
4. The thermometer will switch-off after 5 seconds of idling and retain the selected temperature unit.




FIGURE 6. CHANGING TEMPERATURE SCALES

Cleaning and Storage

This device does not require any special maintenance. If required, use a cotton swab to clean the thermometer's probe if you see contamination on it.

Protect the device from water or other liquids. When not in use, store the thermometer in its original packaging in a cool, dry place and out of the reach of children.

Recharging the Battery

This thermometer is manufactured with a rechargeable battery. If the battery icon (⎓) begins to flash, the battery requires recharging.

Insert the charging cable into the thermometer's micro USB port (see Fig. 7).

Insert the other end of the charging cable into a USB wall socket adaptor and an appropriate power point. Allow a minimum of 2 hours to charge.




FIGURE 7. RECHARGING THE BATTERY

Error Messages

Error Message	Problem	Solution
HI	Device temperature is $\geq 42.2^{\circ}\text{C}$ (108°F)	Ensure the device is being used within the $10^{\circ}\text{C} - 40^{\circ}\text{C}$ temperature operating conditions range noted in the specifications.
LO	Device temperature is $\leq 32^{\circ}\text{C}$ (89.6°F)	Place the device in an environment where temperature is between $10^{\circ}\text{C} - 40^{\circ}\text{C}$ to equilibrate, and wait 15 minutes before using.
ERR	Sensor fault	Contact the Airssential Customer Care Centre on (02) 9708 5560.
ERRE	Internal device fault	
⎓	Low battery	Recharge battery

Warranty

The LifeTemp® Touch Infrared Forehead Thermometer is covered by a limited warranty offered by Airssential Home Health Care Solutions ("Airssential"), a division of Boian Surgical Pty Ltd., to the original purchaser of this product. Proof of purchase is required to substantiate warranty claims. Airssential warrants that for a period of two years from the date of purchase, this product will be free from defects in material and workmanship.

Airssential will replace the device if it is defective during the warranty period, provided the product's use has been in accordance with the instructions provided. If the product is no longer available, replacement may be made with a similar product. Airssential shall not be liable for any incidental or consequential damages caused by the breach of any express or implied warranty.

This warranty excludes replacement due to improper use, accidental user damage or the device's battery which is warranted to be free from defects in workmanship and materials at the time of delivery.

Service

Do not disassemble the device as it cannot be serviced at home. Unauthorised disassembly of this device will void the warranty.

If service is required, contact the Airssential Customer Care Service Centre on (02) 9708 5560. Send this device, together with proof of purchase and a cheque or money order for \$9.90 to cover return shipping and insurance to the Airssential Customer Care Service Centre, 122 Gow Street, Padstow NSW 2211. Please include your name, address and a contact telephone number with your request.

Disposal

Waste medical products should not be disposed of with household rubbish. If this thermometer is no longer required or it is no longer serviceable, dispose of the device and its accessories in compliance with local regulations. Please recycle where facilities exist.

Specifications

Measurement Mode	Forehead Mode
Measurement Unit	Celsius ($^{\circ}\text{C}$) or Fahrenheit ($^{\circ}\text{F}$)
Backlight	Blue Backlight
Operating Conditions	$10-40^{\circ}\text{C}$ ($50-104^{\circ}\text{F}$), RH 10-85%
Storage Conditions	$20-50^{\circ}\text{C}$ ($-4-122^{\circ}\text{F}$), RH 10-85%
Forehead measurement Range and Accuracy	Range: $32 - 42^{\circ}\text{C}$ ($89.6-108.0^{\circ}\text{F}$) Accuracy: $\pm 0.2^{\circ}\text{C}$ ($\pm 0.4^{\circ}\text{F}$)
Other Range	$\pm 0.3^{\circ}\text{C}$ ($\pm 0.5^{\circ}\text{F}$)
Display	LCD
Memory	Last 12 measurements
Auto Power-Off	Automatic switch-off after 60 seconds of idling
Dimensions	51 x 46 x 46 mm
Weight	45g


Distributed By:

Airssential Home Health Care Solutions

122 Gow Street, Padstow NSW 2211

Phone (02) 9708 5560

Airssential® Home Health Care Solutions is a division of Boian Surgical Pty Ltd. Airssential® and LifeTemp® are the registered trademarks of Boian Surgical Pty Ltd. These products are distributed under a quality management system certified as complying with ISO9001:2015.



MADE IN CHINA